

Containers & Cups

● Disposable cups and plastic containers.



300 ml container



100 ml disposable cup



200 ml disposable cup



500 ml disposable cup



150 ml container



250 ml container



300 ml anti-UV container



360 ml container [Type K] With keyways



550 ml container



650 ml container



750 ml container



1000 ml container



4000 ml container



300 ml SUS
300 ml SUS container

*Other types of containers (eg. container lid with a hole) are also available.
Contact us for details.

Correspondence Table

○: Adapter required.

container \ compatible model	AR-100	ARE-310 ARE-250 AR-250	ARE-400TWIN	ARE-500 AR-500	ARV-310 ARV-310LED	ARV-5000 ARV-3000TWIN	SR-500 SP-500	ARV-930TWIN
100ml disposable cup	○	○	○	○	○			
200ml disposable cup		○	○	○	○			
500ml disposable cup				○				○
150 ml container		○	○	○	○	○	○	
250 ml container		●	●	○	●	○	●	
300 ml container		●	●	○	●	○	●	○
300ml anti-UV container		●	●	○	●	○	●	○
360 ml container				○		○		○
550 ml container				●		○		○
650 ml container				●		○		
750 ml container						○		●
1000 ml container						○		
4000 ml container						●		